





Appendix G – Result Interpretation Quick Reference

Result	Interpretation
<p>The control line (C) and the test line (T) are present.</p> 	<p><input checked="" type="checkbox"/> Positive <input type="checkbox"/> Negative <input type="checkbox"/> Invalid</p>
<p>The control line (C) but the test line (T) is faint.</p> 	<p><input checked="" type="checkbox"/> Positive <input type="checkbox"/> Negative <input type="checkbox"/> Invalid</p>
<p>Only the control line (C) is present.</p> 	<p><input type="checkbox"/> Positive <input checked="" type="checkbox"/> Negative <input type="checkbox"/> Invalid</p>
<p>The control line (C) is NOT present</p> 	<p><input type="checkbox"/> Positive <input type="checkbox"/> Negative <input checked="" type="checkbox"/> Invalid</p>